

### **Product datasheet for TA351910**

## **VAMP2 Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type: Primary Antibodies** 

**Applications:** IHC, WB

Recommended Dilution: WB: 500-2000

WB positive control: A172, Jurkat and Hela cells, Mouse brain and Human fetal brain tissue

IHC: 25-100

Positive control: Human thyroid cancer Predicted cell location: Cytoplasm

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: lgG

Clonality: Polyclonal

Immunogen: Synthetic peptide of human VAMP2

pH7.4 PBS, 0.05% NaN3, 40% Glyceroln Formulation:

Concentration: lot specific

**Purification:** Antigen affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stable for 12 months from date of receipt. Stability:

**Predicted Protein Size:** 13 kDa

Gene Name: vesicle associated membrane protein 2

Database Link: NP 055047

Entrez Gene 6844 Human

P63027



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



#### Background:


The protein encoded by this gene is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family. Synaptobrevins/VAMPs, syntaxins, and the 25-kD synaptosomal-associated protein SNAP25 are the main components of a protein complex involved in the docking and/or fusion of synaptic vesicles with the presynaptic membrane. This gene is thought to participate in neurotransmitter release at a step between docking and fusion. The protein forms a stable complex with syntaxin, synaptosomal-associated protein, 25 kD, and synaptotagmin. It also forms a distinct complex with synaptophysin. It is a likely candidate gene for familial infantile myasthenia (FIMG) because of its map location and because it encodes a synaptic vesicle protein of the type that has been implicated in the pathogenesis of FIMG.

Synonyms: SYB2; VAMP-2

**Protein Families:** Druggable Genome, Secreted Protein, Transmembrane

**Protein Pathways:** SNARE interactions in vesicular transport

# **Product images:**



Gel: 12%SDS-PAGE Lysate: 40 µg Lane 1-5: A172 cells Jurkat cells Hela cells

Mouse brain tissue Human fetal brain tissue

Primary antibody: TA351910 (VAMP2 Antibody) at


dilution 1/200

Secondary antibody: Goat anti rabbit IgG at


1/8000 dilution

Exposure time: 30 seconds






Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA351910 (VAMP2 Antibody) at dilution 1/40 (Original magnification: ×200)




Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using TA351910 (VAMP2 Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using TA351910 (VAMP2 Antibody) at dilution 1/40 (Original magnification: ×200)





Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using TA351910 (VAMP2 Antibody) at dilution 1/40, treated with synthetic peptide. (Original magnification: ×200)