

## Product datasheet for **TA352767**

### IL4 Mouse Monoclonal Antibody [Clone ID: DL-8]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	DL-8
Recommended Dilution:	ELISA, ELISPOT, Neutralization assay
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	IL-4
Formulation:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Concentration:	lot specific
Purification:	Ion exchange chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	16.2 kDa
Gene Name:	interleukin 4
Database Link:	<a href="#">NP_958427</a> <a href="#">Entrez Gene 287287 Rat P20096</a>
Background:	The neutralizing mouse anti-rat IL-4 monoclonal antibody (clone DL-8) can be used in ELISA and ELISPOT assays.
Synonyms:	BCGF-1; BCGF1; Binetrakin; BSF-1; BSF1; IL-4; MGC79402; Pitrakinra



[View online »](#)