

Product datasheet for **TA352769**

Ifnb1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	ELISA, Neutralization assay, Western Blot
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	For immunization recombinant rat IFN-beta (CHO-derived) is used
Formulation:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Concentration:	lot specific
Purification:	Antibodies were sequentially purified by ammonium sulphate precipitation and protein A-affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	22.1 kDa
Gene Name:	interferon beta 1
Database Link:	NP_062000 Entrez Gene 24481 Rat P70499
Background:	The neutralizing rabbit anti-rat IFN-beta polyclonal antibody can be used in ELISA and Western Blot. The antibody is affinity purified.
Synonyms:	IFB; IFF; IFN-beta; IFNB; MGC96956



[View online »](#)