

## Product datasheet for **TA352772**

### IL4 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	ELISA, ELISPOT
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	For immunization recombinant rat IL-4 (E.coli-derived) is used
Formulation:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Concentration:	lot specific
Purification:	Antibodies were sequentially purified by ammonium sulphate precipitation and protein A-affinity chromatography.
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	16.2 kDa
Gene Name:	interleukin 4
Database Link:	<a href="#">NP_958427</a> <a href="#">Entrez Gene 287287 Rat</a> <a href="#">P20096</a>
Background:	The biotinylated rabbit anti-rat IL-4 polyclonal antibody can be used in ELISA and ELISPOT assays. The antibody is affinity purified.
Synonyms:	BCGF-1; BCGF1; Binetrakin; BSF-1; BSF1; IL-4; MGC79402; Pitrakinra



[View online »](#)