

## Product datasheet for **TA352777**

### IL2 Rat Monoclonal Antibody [Clone ID: JES6-5H4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	JES6-5H4
Recommended Dilution:	ELISA, Flow cytometry, Immunoprecipitation, Immunohistochemistry, Neutralization assay, Western Blot
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	IL-2
Formulation:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Concentration:	lot specific
Purification:	Protein G affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	19.9 kDa
Gene Name:	interleukin 2
Database Link:	<a href="#">NP_032392</a> <a href="#">Entrez Gene 16183 Mouse</a> <a href="#">P04351</a>
Background:	The neutralizing monoclonal antibody to mouse IL-2 (clone JES6-5H4) has shown to be suitable for ELISA (Carson R.T. et al. 1999 and Gorczyński R.M. et al. 1999), flow cytometry, Western Blot and immunoprecipitation.
Synonyms:	aldesleukin; IL-2; interleukin-2; lymphokine; TCGF



[View online »](#)