

Product datasheet for **TA352778**

IL4 Rat Monoclonal Antibody [Clone ID: 11B11]

Product data:

Product Type:	Primary Antibodies
Clone Name:	11B11
Recommended Dilution:	ELISA, ELISPOT, Intracellular staining, Neutralization assay, Western Blot
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	IL-4
Formulation:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Concentration:	lot specific
Purification:	Ion exchange chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	16.3 kDa
Gene Name:	interleukin 4
Database Link:	NP_067258 Entrez Gene 16189 Mouse P07750
Background:	The monoclonal antibody to mouse IL-4 (clone 11B11) has shown to be suitable for intracellular staining (Assenmacher M. et al. 1998), ELISA (Carson R.T. et al. 1999) and Western Blot. The 11B11 antibody can also successfully be used in immunohistochemistry and neutralization assays (Chaudhri G. et al. 2004 and Hammond K.J. et al. 1998).
Synonyms:	BCGF-1; BCGF1; Binetrakin; BSF-1; BSF1; IL-4; MGC79402; Pitrakinra



[View online »](#)