

Product datasheet for **TA352779**

IL4 Rat Monoclonal Antibody [Clone ID: BVD6-24G2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BVD6-24G2
Recommended Dilution:	ELISA, Flow cytometry
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant mouse IL-4 (E.coli)
Formulation:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Concentration:	lot specific
Purification:	Protein G affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	16.3 kDa
Gene Name:	interleukin 4
Database Link:	NP_067258 Entrez Gene 16189 Mouse P07750
Background:	The monoclonal antibody to mouse IL-4 (clone BVD6-24G2) can be used in ELISA (Chatelain R. et al. 1992) and flow cytometry.
Synonyms:	BCGF-1; BCGF1; Binetrakin; BSF-1; BSF1; IL-4; MGC79402; Pitrakinra



[View online »](#)