

Product datasheet for **TA352780**

IL6 Rat Monoclonal Antibody [Clone ID: MP5-20F3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	MP5-20F3
Recommended Dilution:	ELISA, Intracellular staining, Neutralization assay, Western Blot
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant mouse IL-6
Formulation:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Concentration:	lot specific
Purification:	Protein G chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	24.8 kDa
Gene Name:	interleukin 6
Database Link:	NP_112445 Entrez Gene 16193 Mouse P08505
Background:	The neutralizing rat anti-mouse IL-6 monoclonal antibody (clone MP5-20F3) can be used in ELISA and Western Blot.
Synonyms:	BSF-2; BSF2; CDF; HGF; HSF; IFNB2; IL-6; OTTHUMP00000158544



[View online »](#)