

Product datasheet for **TA352782**

Ifng Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	Neutralization assay, Western Blot
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	For immunization recombinant mouse IFN-gamma (CHO-derived) is used
Formulation:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Concentration:	lot specific
Purification:	Antibodies were sequentially purified by ammonium sulphate precipitation and protein A-affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	18.4 kDa
Gene Name:	interferon gamma
Database Link:	NP_032363 Entrez Gene 15978 Mouse P01580
Background:	The neutralizing polyclonal antibody to mouse IFN-gamma is suitable for Western Blot. The antibody is affinity purified.
Synonyms:	IFG; IFI; IFN-gamma



[View online »](#)