

Product datasheet for **TA352784**

Tnf Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	ELISA, ELISPOT
Reactivity:	Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	For immunization recombinant mouse TNF-alpha (E.coli-derived) is used
Formulation:	Prior to lyophilization: 0.5 ml PBS + 125 mM trehalose.
Concentration:	lot specific
Purification:	Antibodies were sequentially purified by ammonium sulphate precipitation and protein A-affinity chromatography.
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	24.7 kDa
Gene Name:	tumor necrosis factor
Database Link:	NP_001265530 Entrez Gene 24835 Rat Entrez Gene 21926 Mouse P06804
Background:	The biotinylated polyclonal antibody to mouse TNF-alpha is suitable for ELISA and ELISPOT assays. The antibody is affinity purified.
Synonyms:	Cachectin; DIF; TNF-a; TNF-alpha; TNFA; TNFSF2



[View online »](#)