

Product datasheet for **TA353974**

Rabbit Monoclonal Antibody [Clone ID: SP137]

Product data:

Product Type:	Primary Antibodies
Clone Name:	SP137
Recommended Dilution:	IHC-P
Host:	Rabbit
Clonality:	Monoclonal
Formulation:	PBS/1% BSA buffer pH 7.6 with less than 0.1% sodium azide (concentration 200g/mL).
Purification:	protein A/G
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Background:	The Negative Control is a recombinant rabbit monoclonal antibody against a small molecule, which is normally not present in human or animal. The product will not cross-react with human tissues when used at the recommended concentration.



[View online »](#)