

## Product datasheet for **TA354189**

### Chromogranin A (CHGA) Rabbit Polyclonal Antibody

#### Product data:

|                       |   |
|-----------------------|---|
| Product Type:         | Primary Antibodies  |
| Recommended Dilution: | WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml   |
| Reactivity:           | Human, Rat, Monkey  |
| Host:                 | Rabbit  |
| Isotype:              | IgG   |
| Clonality:            | Polyclonal  |
| Formulation:          | This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.  |
| Purification:         | The Rabbit IgG is purified by Epitope Affinity Purification   |
| Conjugation:          | Unconjugated  |
| Storage:              | Store at -20°C as received.   |
| Stability:            | Stable for 12 months from date of receipt.  |
| Gene Name:            | chromogranin A  |
| Database Link:        | <a href="#">NP_001266</a><br><a href="#">Entrez Gene 24258 Rat</a> <a href="#">Entrez Gene 574388 Monkey</a> <a href="#">Entrez Gene 1113 Human</a><br><a href="#">P10645</a> |
| Synonyms:             | CGA   |
| Protein Families:     | Druggable Genome, Secreted Protein  |



[View online »](#)