

Product datasheet for **TA354231**

Parathyroid Hormone (PTH) Rabbit Polyclonal Antibody

Product data:

| | |
|-----------------------|--|
| Product Type: | Primary Antibodies |
| Recommended Dilution: | WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml |
| Reactivity: | Human |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Formulation: | This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer. |
| Purification: | The Rabbit IgG is purified by Epitope Affinity Purification |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Gene Name: | parathyroid hormone |
| Database Link: | NP_000306 Entrez Gene 5741 Human P01270 |
| Synonyms: | PTH1 |
| Protein Families: | Druggable Genome, Secreted Protein |



[View online »](#)