

## Product datasheet for **TA354252**

### Growth Hormone (GH1) Mouse Monoclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
Reactivity:	Human
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Formulation:	This purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The Rabbit IgG is purified by Protein A-Affinity Chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	growth hormone 1
Database Link:	<a href="#">NP_000506</a> <a href="#">Entrez Gene 2688 Human</a> <a href="#">P01241</a>
Synonyms:	GH; GH-N; GHN; hGH-N; IGHD1B
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway, Neuroactive ligand-receptor interaction



[View online »](#)