

## Product datasheet for **TA354258**

### Bak (BAK1) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
Reactivity:	Human, Mouse, Hamster
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The Rabbit IgG is purified by Epitope Affinity Purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	BCL2 antagonist/killer 1
Database Link:	<a href="#">NP_001179</a> <a href="#">Entrez Gene 12018 Mouse</a> <a href="#">Entrez Gene 578 Human</a> <a href="#">Q16611</a>
Synonyms:	BAK; BAK-LIKE; BCL2L7; CDN1
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Transmembrane



[View online »](#)