

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA354393

HSF1 Rat Monoclonal Antibody [Clone ID: SPM401]


Product data:

Product Type:	Primary Antibodies
Clone Name:	SPM401
Applications:	IHC, WB
Recommended Dilution:	WB 0.1-1 μg/ml ELISA 0.01-0.1 μg/ml IP 2-5 μg/ml IHC 2-10 μg/ml FC 5-10 μg/ml
Reactivity:	Human, Mouse, Rat
Host:	Rat
lsotype:	lgG1
Clonality:	Monoclonal
Formulation:	This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The mouse IgG is purified by Protein A-Affinity Chromatography according to Isotyping
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	heat shock transcription factor 1
Database Link:	<u>NP_005517</u> <u>Entrez Gene 15499 MouseEntrez Gene 79245 RatEntrez Gene 3297 Human</u> <u>Q00613</u>
Synonyms:	HSTF1
Protein Families:	Transcription Factors




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



WB: The cell lysate was derived from MCF7 was resolved onto 7.5% SDSPAGE, transferred onto NC membrane, and immunoprobed by Rat anti-HSF1 antibody (at 1:500 dilution). A major immunoreactive band is observed around 70kDa.



IHC: Human breast carcinoma tissue stained with Rat anti-HSF1 antibody, at 1:100 for 10 min at RT. Staining of formalin-fixed tissue requires boiling tissue sections in 10 mM Citrate Buffer, pH 6.0 for 10 min followed by cooling at RT for 20 min.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US