

Product datasheet for **TA354449**

IRAKM (IRAK3) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The Rabbit IgG is purified by Epitope Affinity Purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	interleukin 1 receptor associated kinase 3
Database Link:	NP_001135995 Entrez Gene 73914 Mouse Entrez Gene 314870 Rat Entrez Gene 11213 Human Q9Y616
Synonyms:	ASRT5; IRAKM
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Apoptosis, Neurotrophin signaling pathway



[View online »](#)