

Product datasheet for **TA354546**

ADAM8 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	A synthetic peptide corresponding to the human CD156b at Pro- domain. This sequence is identical among human, mouse and/or rat origins.
Formulation:	This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The Rabbit IgG is purified by Epitope Affinity Purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	~100 kDa
Gene Name:	ADAM metallopeptidase domain 8
Database Link:	NP_001100 Entrez Gene 11501 Mouse Entrez Gene 101 Human P78325



[View online »](#)

Background:

CD156b is a type I transmembrane glycoprotein belong to the ADAM (a disintegrin and metalloprotease domain) family. This single-pass type-1 824 aa transmembrane glycoprotein contains a signal sequence, an extracellular domain (pro-domain, catalytic domain, disintegrin-cysteine-rich domain, glycosylation sites), transmembrane and cytoplasmic domain. CD156b is a multi-domain protein that includes a Zn-dependent protease domain in the extracellular portion. It releases the soluble forms of TNF and TGF α from cells. CD156b is identified by its ability to cleave the transmembrane form of TNF between Ala-76 and Val-77 to generate the soluble form. CD156b is therefore a prospective therapeutic target in human cancer. CD156b is expressed on T cells, neutrophils, endothelial cells, monocytes, dendritic cells, macrophages, polymorphonuclear leukocytes and myocytes. And non-lymphoid tissues including heart, placenta, brain, muscle.

Synonyms:

CD156; CD156a; MS2

Protein Families:

Druggable Genome, Transmembrane

Product images:

WB: The cell lysate derived from Jurkat was immunoprobed by Rabbit anti-CD156b (NT) at 1:500. An immunoreactive band is observed around ~