

Product datasheet for TA354599

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SOX4 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Recommended Dilution: WB 0.1-1 μg/ml ELISA 0.01-0.1 μg/ml IP 2-5 μg/ml IHC 2-10 μg/ml FC 5-10 μg/ml

Reactivity: Human
Host: Rabbit
Isotype: IgG

Clonality: Polyclonal

Immunogen: A synthetic peptide corresponding to the c-terminus (from 720AA-755aa) of human sox6

protein. This sequence is identical in human, mouse, rat, bovine, canis and chicken species.

Formulation: This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2)

containing antibody stabilizer.

Purification: The Rabbit IgG is purified by Epitope Affinity Purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: SRY-box 4

Database Link: NP 003098

Entrez Gene 6659 Human

Q06945

Background: SOX6 gene encodes a member of the D subfamily of sex determining region y-related

transcription factors that are characterized by a conserved DNA-binding domain termed the high mobility group box and by their ability to bind the minor groove of DNA. SOX6 protein is a transcriptional activator that is specifically expressed in the developing central nervous system and in the early stages of chondrogenesis in mouse embryos. It is found to be essential to both cartilage formation/chondrogenesis and obesity-related insulin resistance.

The SOX6 protein is also detected in glioma patient sera and brain tissues.

Synonyms: EVI16

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Stem cell

relevant signaling - TGFb/BMP signaling pathway, Transcription Factors

