

## Product datasheet for **TA354632**

### Annexin A1 (ANXA1) Rabbit Polyclonal Antibody

#### Product data:

|                       |  |
|-----------------------|--|
| Product Type:         | Primary Antibodies   |
| Recommended Dilution: | WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml  |
| Reactivity:           | Human  |
| Host:                 | Rabbit   |
| Isotype:              | IgG  |
| Clonality:            | Polyclonal   |
| Formulation:          | This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer. |
| Purification:         | The Rabbit IgG is purified by Epitope Affinity Purification  |
| Conjugation:          | Unconjugated   |
| Storage:              | Store at -20°C as received.  |
| Stability:            | Stable for 12 months from date of receipt.   |
| Gene Name:            | annexin A1   |
| Database Link:        | <a href="#">NP_000691</a><br><a href="#">Entrez Gene 301 Human</a><br><a href="#">P04083</a>                             |
| Synonyms:             | ANX1; LPC1   |
| Protein Families:     | Druggable Genome   |



[View online »](#)