

Product datasheet for **TA354636**

Calca Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
Reactivity:	Human, Mouse, Rat, Monkey, Dog
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The Rabbit IgG is purified by Epitope Affinity Purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	calcitonin/calcitonin-related polypeptide, alpha
Database Link:	NP_001029126 Entrez Gene 796 Human Entrez Gene 24241 Rat Entrez Gene 403946 Dog Entrez Gene 700649 Monkey Entrez Gene 12310 Mouse Q99JA0
Synonyms:	CALC1; CGRP; CGRP-I; CGRP1; CT; katalcalcin; KC; MGC126648; PCT; Procalcitonin



[View online »](#)