

## Product datasheet for **TA354762**

### alpha 1 Fetoprotein (AFP) Mouse Monoclonal Antibody [Clone ID: ABM001]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	ABM001
Recommended Dilution:	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The mouse IgG is purified by Protein A-Affinity Chromatography according to Isotyping
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	alpha fetoprotein
Database Link:	<a href="#">NP_001125</a> <a href="#">Entrez Gene 174 Human</a> <a href="#">P02771</a>
Synonyms:	AFPD; FETA; HPAFP
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)