

## Product datasheet for **TA354794**

### Lhb Mouse Monoclonal Antibody [Clone ID: ABM501]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	ABM501
Recommended Dilution:	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The mouse IgG is purified by Protein A-Affinity Chromatography according to Isotyping
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	luteinizing hormone beta
Database Link:	<a href="#">NP_032523</a> <a href="#">Entrez Gene 3972 Human</a>
Synonyms:	CGB4; hLHB; LH-B; LSH-B; LSH-beta



[View online »](#)