

## Product datasheet for **TA354801**

### **tat Mouse Monoclonal Antibody**

#### **Product data:**

|                              |  |
|------------------------------|--|
| <b>Product Type:</b>         | Primary Antibodies   |
| <b>Recommended Dilution:</b> | WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml  |
| <b>Reactivity:</b>           | Human  |
| <b>Host:</b>                 | Mouse  |
| <b>Isotype:</b>              | IgG1/IgM   |
| <b>Clonality:</b>            | Monoclonal   |
| <b>Immunogen:</b>            | Tat rec. prot.   |
| <b>Formulation:</b>          | This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer. |
| <b>Purification:</b>         | The mouse IgG is purified by Protein A-Affinity Chromatography according to Isotyping                                    |
| <b>Conjugation:</b>          | Unconjugated   |
| <b>Storage:</b>              | Store at -20°C as received.  |
| <b>Stability:</b>            | Stable for 12 months from date of receipt.   |
| <b>Database Link:</b>        | <a href="#">Entrez Gene 155871</a>   |



[View online »](#)