

## Product datasheet for **TA354892**

### **RALA Rabbit Polyclonal Antibody**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Recommended Dilution:</b>	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Formulation:</b>	This affinity purified antibody is supplied in sterile Tris-buffered saline (pH7.2) containing antibody stabilizer.
<b>Purification:</b>	The Rabbit IgG is purified by Epitope Affinity Purification
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	RALA Ras like proto-oncogene A
<b>Database Link:</b>	<a href="#">NP_005393</a> <a href="#">Entrez Gene 56044 Mouse</a> <a href="#">Entrez Gene 81757 Rat</a> <a href="#">Entrez Gene 5898 Human</a> <a href="#">P11233</a>
<b>Synonyms:</b>	RAL
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Pancreatic cancer, Pathways in cancer



[View online »](#)