

Product datasheet for **TA354912**

AXL Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	This affinity purified antibody is supplied in sterile Tris-buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The Rabbit IgG is purified by site-modified Epitope Affinity Purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	AXL receptor tyrosine kinase
Database Link:	NP_001265528 Entrez Gene 26362 Mouse Entrez Gene 308444 Rat Entrez Gene 558 Human P30530
Synonyms:	ARK; JTK11; Tyro7; UFO
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane



[View online »](#)