

Product datasheet for **TA354918**

TJP2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Recommended Dilution:	WB 0.1-1 µg/ml ELISA 0.01-0.1 µg/ml IP 2-5 µg/ml IHC 2-10 µg/ml FC 5-10 µg/ml
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The Rabbit IgG is purified by Epitope Affinity Purification
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	tight junction protein 2
Database Link:	NP_001163885 Entrez Gene 21873 MouseEntrez Gene 115769 RatEntrez Gene 9414 Human Q9UDY2
Synonyms:	C9DUPq21.11; DFNA51; DUP9q21.11; PFIC4; X104; ZO2
Protein Pathways:	Tight junction, Vibrio cholerae infection



[View online »](#)