

## Product datasheet for **TA354950**

### **TOR2A Rabbit Polyclonal Antibody**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Formulation:</b>	This affinity purified antibody is supplied in sterile Tris-buffered saline (pH7.2) containing antibody stabilizer.
<b>Purification:</b>	The Rabbit IgG is purified by site-modified Epitope Affinity Purification.
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	torsin family 2 member A
<b>Database Link:</b>	<a href="#">NP_001078816</a> <a href="#">Entrez Gene 30933 Mouse</a> <a href="#">Entrez Gene 362112 Rat</a> <a href="#">Entrez Gene 27433 Human</a> <a href="#">Q5JU69</a>
<b>Synonyms:</b>	TORP1
<b>Protein Families:</b>	Secreted Protein, Transmembrane



[View online »](#)