

Product datasheet for **TA354952**

CNPY4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	This affinity purified antibody is supplied in sterile Tris-buffered saline (pH7.2) containing antibody stabilizer.
Purification:	The Rabbit IgG is purified by site-modified Epitope Affinity Purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	canopy FGF signaling regulator 4
Database Link:	NP_689968 Entrez Gene 66455 Mouse Entrez Gene 363886 Rat Entrez Gene 245812 Human Q8N129
Synonyms:	PRAT4B



[View online »](#)