

Product datasheet for TA355071

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PD-L1 (CD274) Mouse Monoclonal Antibody [Clone ID: 1D7]

Product data:

Product Type: Primary Antibodies

Clone Name: 1D7

Applications: IF, IHC, WB

Recommended Dilution: WB: 0.5-1µg/mL.IHC starting at 2-5µg/mL.IF start at 5µg/mL.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: PD-L1 antibody was raised against the extracellular domain of human PD-L1.

Formulation: PD-L1 Antibody is supplied in PBS containing 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: PD-L1 Antibody is supplied as protein A purified IgG1.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: Predicted: 32 kDa; Observed: 45 kDa

Gene Name: CD274 molecule

Database Link: NP 054862

Entrez Gene 29126 Human

Q9NZQ7



Background:

PD-L1 Antibody: Cell-mediated immune responses are initiated by T lymphocytes that are themselves stimulated by cognate peptides bound to MHC molecules on antigen-presenting cells (APC) (1). T-cell activation is generally self-limited as activated T cells express receptors such as PD-1 (also known as PDCD-1) that mediate inhibitory signals from the APC (2). PD-1 can bind two different but related ligands, PD-L1 and PD-L2. PD-L1 is a B7-related protein that inhibits cell-mediated immune responses by reducing the secretion of IL-2 and IL-10 from memory T cells (3). This suggests that PD-L1 may be useful in reducing allogenic CD4+ memory T-cell responses to endothelial cells, thereby reducing the likelihood of host immune responses to allografts. PD-L1 also functions as an immune checkpoint protein, and multiple anti-PD-L1 antibodies are currently in phase II and III clinical trials, with one antibody already approved for the treatment of cancer (4).


Synonyms:

B7-H; B7-H1; B7H1; MGC142294; MGC142296; PD-L1; PDCD1L1; PDCD1LG1; PDL1

Note:


PD-L1 antibody can be used for detection of PD-L1 by Western blot at 0.5 - 1 μ g/mL. Antibody can also be used for immunohistochemistry starting at 2 - 5 μ g/mL. For immunofluorescence start at 5 μ g/mL.

Product images:




Western blot analysis of PD-L1 in overexpressing HEK293 cells PD-L1 antibody at 0.25 and 0.5ug/ml






Immunohistochemistry of PD-L1 in human stomach carcinoma tissue with PD-L1 antibody at 5ug/ml.




Immunohistochemistry of PD-L1 in human tonsil tissue with PD-L1 antibody at 5ug/ml.




Immunofluorescence of PD-L1 in transfected HEK293 cells with PD-L1 antibody at 2ug/ml.






Immunofluorescence of PD-L1 in human stomach carcinoma tissue with PD-L1 antibody at 2ug/ml.



Immunofluorescence of PD-L1 in human tonsil tissue with PD-L1 antibody at 2ug/ml.



Immunocytochemistry of PD-L1 in transfected HEK293 cells with PD-L1 antibody at 1ug/ml. Lower left: Immunocytochemistry in transfected HEK293 cells with control mouse IgG antibody at 1ug/ml.