

Product datasheet for **TA355150**

SARS Matrix protein antibody Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	SARS-CoV
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Raised against a synthetic peptide corresponding to amino acids at the amino-terminus of the SARS matrix protein.
Formulation:	PBS containing 0.02% sodium azide
Concentration:	lot specific
Purification:	Affinity chromatography purified via peptide column
Conjugation:	Unconjugated
Storage:	Can be stored at -20 °C. Avoid repeated freeze thaw cycles.
Stability:	one year from dispatch
Database Link:	P59596
Background:	A novel coronavirus has recently been identified as the causative agent of SARS (Severe Acute Respiratory Syndrome). Coronaviruses are a major cause of upper respiratory diseases in humans. The genomes of these viruses are positive-stranded RNA approximately 27-31kb in length. The M protein (Membrane protein, Matrix protein) is one of the major structural viral proteins. It is an integral membrane protein involved in the budding of the viral particles and interacts with S (Spike) protein and the nucleocapsid protein.



[View online »](#)