

Product datasheet for **TA355159**

SARS-CoV-2 Spike Rabbit Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, FC, IF, IHC, IP
Recommended Dilution:	Elisa: 1:1000 - 1:5000 FCM: 1:50 - 1:200 IHC: 1:50 - 1:200 IF: 1:50 - 1:200 IP: 1:50 - 1:200
Reactivity:	SARS-CoV-2
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	recombinant protein
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	1 year
Gene Name:	S Protein
Database Link:	Entrez Gene 43740568 SARS-CoV-2



[View online »](#)