

Product datasheet for TA355211

Human IgA (Secretory component, free and bound) Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA (Secretory component, free and bound) Sheep Polyclonal Antibody
Applications:	E, ICC, IHC, WB
Recommended Dilution:	This immunoconjugate is not pre-diluted. The optimum working dilution of each conjugate should be established by titration before being used. Excess labelled antibody must be avoided because it may cause high unspecific background staining and interfere with the specific signal.
Reactivity:	Human
Host:	Sheep
Immunogen:	Purified free human secretory component isolated from pooled milk is used for immunization. Freund's complete adjuvant is used in the first step of the immunization procedure.
Isotype:	IgG
Formulation:	Horseradish peroxidase-coupled purified hyperimmune sheep IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.
Reconstitution Method:	The antibody is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divided into small aliquots, frozen and stored at or below -20°C.
Conjugation:	HRP
Storage:	The lyophilized conjugate can be stored at +4°C short term; prolonged storage at or below -20°C.



[View online »](#)