

Product datasheet for **TA355217**

alpha Sarcoglycan (SGCA) Mouse Monoclonal Antibody [Clone ID: Ad1/20A6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	Ad1/20A6
Applications:	IHC
Recommended Dilution:	1:100–1:200
Reactivity:	Human, Hamster, Mouse, Pig, Rabbit, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Fusion protein containing amino acids 217 to 289 of the rabbit adhalin sequence
Specificity:	Human alpha-sarcoglycan, also known as adhalin. Also crossreacts strongly with alpha-sarcoglycan in sections of muscle from mouse, rat, rabbit, hamster and pig. Does not react with chicken muscle.
Formulation:	Liquid tissue culture supernatant containing sodium azide as a preservative
Conjugation:	Unconjugated
Storage:	Store at 2-8°C
Stability:	12 months
Gene Name:	sarcoglycan alpha
Database Link:	Entrez Gene 6442 Human Q16586
Background:	Human alpha-sarcoglycan, also known as adhalin. Also crossreacts strongly with alpha-sarcoglycan in sections of muscle from mouse, rat, rabbit, hamster and pig. Does not react with chicken muscle.
Synonyms:	50-DAG; 50DAG; A2; Adhalin; ADL; alpha-sarcoglycan; Alpha-SG; DAG2; DMDA2; Dystroglycan-2; LGMD2D; SCARMD1



[View online »](#)