

Product datasheet for **TA355227**

beta Sarcoglycan (SGCB) Mouse Monoclonal Antibody [Clone ID: **βSarc1/5B1**]

Product data:

Product Type:	Primary Antibodies
Clone Name:	βSarc1/5B1
Applications:	IHC
Recommended Dilution:	1:50–1:200
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Fusion protein RBSG-NT of the human beta-sarcoglycan sequence.
Specificity:	Human beta-sarcoglycan (43 kD)
Formulation:	Liquid tissue culture supernatant containing sodium azide as a preservative
Conjugation:	Unconjugated
Storage:	Store at 2-8°C
Stability:	12 months
Gene Name:	sarcoglycan beta
Database Link:	Entrez Gene 6443 Human Q16585
Background:	Human beta-sarcoglycan (43 kD)
Synonyms:	43DAG; A3b; Beta-SG; LGMD2E; SGC



[View online »](#)