

Product datasheet for **TA355261**

CD117 Mouse Monoclonal Antibody [Clone ID: T595]

Product data:

Product Type:	Primary Antibodies
Clone Name:	T595
Applications:	IHC
Recommended Dilution:	1:40
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Prokaryotic recombinant protein corresponding to the three N-terminal C2-like extracellular domains of the human c-kit oncoprotein molecule
Specificity:	Human c-kit oncoprotein
Formulation:	Liquid tissue culture supernatant containing sodium azide as a preservative
Conjugation:	Unconjugated
Storage:	Store at 2-8°C
Stability:	12 months
Gene Name:	KIT proto-oncogene receptor tyrosine kinase
Database Link:	Entrez Gene 3815 Human P10721
Synonyms:	C-Kit; CD117; PBT; SCFR



[View online »](#)