

Product datasheet for **TA355284**

Dysferlin (DYSF) Mouse Monoclonal Antibody [Clone ID: Ham1/7B6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	Ham1/7B6
Applications:	IHC
Recommended Dilution:	1:20 - 1:40
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Synthetic peptide containing amino acids 1999-2016 of the human dysferlin molecule
Specificity:	Reactive with the dysferlin molecule in human skeletal muscle. Also present in many non-muscle tissues
Formulation:	Lyophilized tissue culture supernatant containing sodium azide as a preservative.
Reconstitution Method:	The user is required to reconstitute the contents of the vial with the correct volume of sterile distilled water as indicated on the vial label
Conjugation:	Unconjugated
Storage:	Store at 2-8°C
Stability:	12 months
Gene Name:	dysferlin
Database Link:	Entrez Gene 8291 Human O75923
Background:	Dysferlin is the protein product of the 2p13 gene that is defective in patients with Limb-Girdle Muscular Dystrophy type 2B (LGMD2B) and Miyoshi Myopathy (MM). Dysferlin is normally localized to the muscle plasma membrane. In patients with LGMD2B and MM, immunoreactivity to dysferlin is severely reduced or lost. Patients with other neuromuscular conditions demonstrate normal labeling patterns.
Synonyms:	dysferlin; FER1L1; FLJ00175; FLJ90168; LGMD2B; OTTHUMP00000202233; OTTHUMP00000202234; OTTHUMP00000202235; OTTHUMP00000202236; OTTHUMP00000202237; OTTHUMP00000202240



[View online »](#)