

Product datasheet for **TA355364**

B7H3 (CD276) Mouse Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Host:	Mouse
Clonality:	Monoclonal
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD276 molecule
Database Link:	Q5ZPR3
Synonyms:	4Ig-B7-H3; B7-H3; B7H3



[View online »](#)