

Product datasheet for **TA355366**

GFP Mouse Monoclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Specificity:	Expected reactivity: Aequorea victoria, Aequorea aequorea, Aequorea forskalea
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	27kDa
Database Link:	AAC53684 P42212
Background:	GFP was originally isolated from the jellyfish



[View online »](#)

Product images:



WB Suggested Anti-GFP Antibody Titration:
1.0ug/ml
Positive Control: rGFP-GST