

## Product datasheet for **TA355410**

### CATSPER3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of human CATSPER3
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	46 kDa
Gene Name:	cation channel sperm associated 3
Database Link:	<a href="#">NP_821138.1</a> <a href="#">Entrez Gene 347732 Human</a> <a href="#">Q86XQ3</a>
Synonyms:	CACRC; MGC126741
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane



[View online »](#)