

Product datasheet for **TA355416**

SLC43A3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human SLC43A3
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	54 kDa
Gene Name:	solute carrier family 43 member 3
Database Link:	NP_001265130.1 Entrez Gene 29015 Human Q8NBI5
Synonyms:	DKFZp762A227; Eeg1; FOAP-13; PRO1659; SEEEG-1
Protein Families:	Transmembrane



[View online »](#)

Product images:

Host: Rabbit
Target Name: SLC43A3
Sample Type: Jurkat Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: SLC43A3
Sample Tissue: Jurkat Whole Cell lysates
Antibody Dilution: 1 μ g/ml