

Product datasheet for **TA355418**

ZBTB8OS Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the Middle region of Human ZBTB8OS
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	20 kDa
Gene Name:	zinc finger and BTB domain containing 8 opposite strand
Database Link:	Entrez Gene 339487 Human Q8IWT0
Background:	The function of this protein remains unknown. There are 2 isoforms produced by alternative splicing.
Synonyms:	ARCH; ARCH2; MGC62007



[View online »](#)

Product images:



Host: Rabbit
Target Name: ZBTB8OS
Sample Tissue: Stomach Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: ZBTB8OS
Sample Type: Stomach Tumor
Antibody Dilution: 1.0ug/ml