

Product datasheet for **TA355428**


FANK1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C-terminal region of Human FANK1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	38 kDa
Gene Name:	fibronectin type III and ankyrin repeat domains 1
Database Link:	NP_660278.3 Entrez Gene 92565 Human Q8TC84
Background:	The function of FANK1 remains unknown.
Synonyms:	1700007B22Rik; HSD13



[View online »](#)

Product images:


Host: Rabbit
Target Name: FANK1
Sample Tissue: 293T Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: FANK1
Sample Type: 293T Whole Cell
Antibody Dilution: 1.0ug/ml