

## Product datasheet for **TA355536**

### Eya2 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Horse, Human, Mouse, Pig, Rabbit, Rat <b>Homology:</b> Cow: 93%; Dog: 79%; Horse: 86%; Human: 100%; Mouse: 93%; Pig: 93%; Rabbit: 93%; Rat: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	58kDa
Gene Name:	EYA transcriptional coactivator and phosphatase 2
Database Link:	<a href="#">NP_569111</a> <a href="#">Entrez Gene 156826 Rat</a> <a href="#">Q6UN47</a>
Background:	Eya2 is a member of the eyes absent (EYA) family of proteins; may play a role in eye development; mouse homolog may act as a transcriptional activator.
Synonyms:	EAB1; MGC10614



[View online »](#)

Product images:



WB Suggested Anti-Eya2 Antibody  
Titration: 1.0 ug/ml  
Positive Control: Rat Muscle