

Product datasheet for **TA355542**


ZNF776 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of Human ZNF776
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	56kDa
Gene Name:	zinc finger protein 776
Database Link:	NP_775903 Entrez Gene 284309 Human Q68DI1
Background:	ZNF776 may be involved in transcriptional regulation.
Synonyms:	DKFZp686K10134; DKFZp781G1213; FLJ38288



[View online »](#)

Product images:


Host: Rabbit
 Target Name: ZNF776
 Sample Tissue: MCF7 Cell Lysate
 Antibody Dilution: 1.0µg/ml

Host: Rabbit
 Target Name: ZNF776
 Sample Type: MCF7 Whole Cell lysates
 Antibody Dilution: 1.0ug/ml