

Product datasheet for **TA355668**

Cnot3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of RAT Cnot3
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Zebrafish Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	57kDa
Gene Name:	CCR4-NOT transcription complex, subunit 3
Database Link:	Entrez Gene 308311 Rat D3ZUV9
Background:	The function of this protein remains unknown.
Synonyms:	KIAA0691; LENG2; NOT3; NOT3H; OTTHUMP00000197713



[View online »](#)

Product images:



Host: Rabbit
Target Name: Cnot3
Sample Tissue: Rat Thymus Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: Cnot3
Sample Type: Rat Thymus lysates
Antibody Dilution: 1.0ug/ml