

Product datasheet for **TA355708**


ZFP2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Guinea Pig, Human Homology: Guinea Pig: 79%; Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	53kDa
Gene Name:	ZFP2 zinc finger protein
Database Link:	NP_085116 Entrez Gene 80108 Human Q6ZN57
Background:	ZFP2 is a probable transcription factor involved in neuronal differentiation and/or phenotypic maintenance.
Synonyms:	FLJ21628; Zfp-2; ZNF751



[View online »](#)

Product images:

WB Suggested Anti-ZFP2 Antibody
Titration: 1.0 ug/ml
Positive Control: Hela Whole Cell