

Product datasheet for **TA355802**

ZNF12 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human ZNF12
Specificity:	Expected reactivity: Human Homology: Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	58kDa
Gene Name:	zinc finger protein 12
Database Link:	Entrez Gene 7559 Human P17014
Background:	ZNF12 is a candidate transcription factor
Synonyms:	GIOT-3; GIOT3; HZF11; KOX3; OTTHUMP00000200756; ZNF325




[View online »](#)

Product images:



WB Suggested Anti-ZNF12 Antibody Titration:
1.25ug/ml
ELISA Titer: 1:1562500
Positive Control: HepG2 cell lysate



Human Intestine