

Product datasheet for **TA355953**


Nkx1-1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of Rat Nkx1-1
Specificity:	Expected reactivity: Mouse, Rat Homology: Mouse: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	40kDa
Gene Name:	NK1 homeobox 1
Database Link:	XP_001065657 Entrez Gene 298981 Rat
Background:	The function of this protein remains unknown.



[View online »](#)

Product images:

Host: Rabbit
Target Name: Nkx1-1
Sample Type: Rat Testis lysates
Antibody Dilution: 1.0ug/ml