

## Product datasheet for **TA355969**

### Aff3 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Horse, Human, Mouse, Pig, Rabbit, Rat, Zebrafish <b>Homology:</b> Cow: 93%; Dog: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Pig: 100%; Rabbit: 100%; Rat: 100%; Zebrafish: 85%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	62kDa
Gene Name:	AF4/FMR2 family, member 3
Database Link:	<a href="#">NP_034808</a> <a href="#">Entrez Gene 16764 Mouse P51827</a>
Background:	Aff3 is a putative transcription activator that may function in lymphoid development and oncogenesis.
Synonyms:	LAF4; MLLT2-like; OTTHUMP00000202880; OTTHUMP00000202881



[View online »](#)

## Product images:



WB Suggested Anti-Aff3 Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Thymus